










# **CAT**<sup>TM</sup> ADJUSTING TOOL

## Indications, Contraindications and Warnings for Use

### Indications for Use


The Chiropractic Adjusting Tool (CAT) is indicated for use when chiropractic and other joint manipulations are appropriate for treatment of pain on any patient. The High Velocity, Low Amplitude (HVLA) adjustment with the CAT reduces muscle resistance providing a more effective treatment than by hand and allows for more control.

### Warnings and Contraindications

-  This device should only be used by medical professionals
-  Do not use on or near the eyes
-  Do not use on or near fractures
-  Do not use on or near open wounds
-  Do not use on or near burned tissue
-  Do not use on patient suffering from severe osteoporosis
-  Do not use for any purpose other than adjustment of skeletal joints or structures
-  Do not use without soft tip
-  If device is damaged or does not perform to specifications; contain device and related materials to avoid continued use and contact JTECH Medical

### **WARNING**

CAT EX is a high thrust device exceeding the thrust force of the CAT. It is intended for use by licensed physicians trained in mechanical thrust instrument use and has a maximum rated output of 52 lb<sub>f</sub>.

  
JTECH Medical  
7633 S Main Building D  
Midvale, UT 84047 USA  
Phone: (385) 695-5000  
www.jtechmedical.com

**EC REP**  
Advena Ltd.  
Tower Business Centre,  
2nd Floor,  
Tower Street, Swatar,  
BKR 4013 Malta

**UK RP**  
Advena Ltd.  
Pure Offices  
Plato Close  
Warwick CV34 6WE  
UK



# **CAT**<sup>TM</sup> ADJUSTING TOOL

## Instructions

### Instructions for Use


1. Place fingers on underside of the lower handle with top handle of CAT against palm of hand insuring the grip on the device is secure.
2. Adjust device by turning the adjusting sleeve. (Minimum force setting is achieved when the adjusting sleeve is closest to the tip. Maximum force setting is achieved when the adjusting sleeve is set to first indicator closest to the sleeve but not tightened against the housing.)
3. Test thrust against your hand and adjust until the desired thrust is set.
4. Ensure CAT fully resets after firing by releasing the finger pads to starting position.



### Instructions for Attaching/Removing Cervical Tip

With a rocking motion, push rubber tip into the base of the cervical tip until it is securely attached yet not difficult to take off.

To remove the cervical tip, pull on the tip while gently rocking it back and forth until it releases.



### Cleaning/Disinfecting Instructions

Wipe the CAT with a soft, dry cloth or CaviWipes<sup>TM</sup> disinfecting towelette. Use CaviWipes only as directed, following all instructions and warnings on packaging. Allow device to dry before use.

 Do not submerge or autoclave.

### Storage

Always store your CAT in a clean, dry area.