



JTECH MEDICAL
MEASURING HUMAN PERFORMANCE



Instrument Registration and Interface Software

QUICK START GUIDE

IMPORTANT! Keep the software in a safe place,
SEPARATE from your computer, in case of computer crash or loss.

Computer System Specifications

For system requirements please visit: <https://www.jtechmedical.com/system-requirements>

Installing IRIS Software

Follow these steps to install your software:

- 1) Insert the JTECH Medical IRIS Software USB Drive into one of your computer's available USB ports.
- 2) Browse the files on the USB drive and double-click on the "IrisSetup.exe" file to launch installer.
- 3) Follow the on-screen instructions to complete the IRIS installation. If you are prompted by User Account Controls, select the "Allow" option.



Registration and License Key

Before using your software you must call JTECH Customer Service to register your software and receive a license key. You will perform the following steps while talking to your Customer Service Representative.

1) Open the IRIS Software

Click on the IRIS icon on your desktop or open the software through the Windows “All Programs” menu.

2) Locate the Site ID and Site Code

The registration screen should automatically appear during installation.

3) Get License Key

There are two ways to obtain a license:

- 1) Click “Request License Via Email” button on registration screen. This will launch your default email program and you should see an email form. Fill in the required additional information and click send. Shortly after, you should receive a license key via email or a phone call.
- 2) Call Customer Service: (385) 695-5000 or (800) 985-8324. Provide the Site ID and Site Code to a representative. They will provide the License Key for you to enter in the field under the License Information tab.

4) Enter the License Key

Type the numbers and letters exactly as provided by the Customer Service Representative into the License Key field.

5) Click “Verify License Key”

6) Click “Continue”

Your registration is complete.

JTECH Software Licensing

JTECH MEDICAL
MEASURING HUMAN PERFORMANCE

If you need to register your software or if you encounter any errors, please contact JTECH Technical Support at 1-800-985-8324 or email us at license@jtechmedical.com.

Product Version

Site ID Site Code

License Information | Advanced | Log

License Type License Quantity

License Issued Options Enabled

License Expires

License Key



License Agreement

As part of the install you are asked to read the IRIS Software license agreement, which contains the terms of use for IRIS Software. Clicking "I Accept" and completing the IRIS Software installation indicates you agree to abide by the licensing agreement. If you do not accept the license agreement, do not click "I Accept" or complete the installation. If you do not accept the license agreement, contact JTECH Customer Service.

Limited Warranty

This software is sold AS IS without further warranty as to performance, merchantability, or fitness for a particular purpose. The entire risk as to the results and performance of this software is assumed by the customer. JTECH warrants the USB on which the software is recorded to be free from defects in materials and workmanship under normal use for a period of 90 days from the date of purchase. If the USB should become defective during the 90-day period, it may be returned to JTECH for replacement without charge. In addition, JTECH warrants the software to be free from significant errors that make it unusable for a period of 90 days from the date of purchase. In the event such errors are found, JTECH will reasonably attempt to correct them, or help the customer avoid them, with efforts JTECH believes suitable to the problem or, at JTECH's option, authorize a refund of the license fee upon return of the program. The customer's sole and exclusive remedy in the event of a defect is expressly limited to the above provision. JTECH makes no warranty against material that has been lost, stolen, or damaged by accident, misused, modified without authorization, or for conflicts that arise out of the use of software not supplied by JTECH. JTECH makes no other warranty express or implied to the customer or any other entity or persons, specifically, JTECH makes no warranty that the software is fit for a particular purpose, any implied warranty of merchantability is expressly and specifically disclaimed.

JTECH will not be liable for special, incidental, consequential, indirect, or other similar damages even if JTECH or its agent has been advised of the possibility of such damages. JTECH is not responsible or liable for damages or costs incurred as a result of loss of time, loss or damage of data, loss of profits or revenue, loss of use of the software, or any other losses whatsoever. In addition, JTECH is not responsible or liable for damages or costs incurred in connection with obtaining substitute software, claims by others, inconvenience or similar costs. In no event will JTECH's liability for any damages to the customer or any other entity or person ever exceed the price paid for the license to use the software, regardless of any form of the claim.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state. The customer and JTECH acknowledge that the software is not intended as a consumer good under state and federal warranty laws.

Additional Information

For more detailed information about IRIS software, visit: <http://manuals.jtechmedical.com#IRIS> to download the IRIS Manual.



JTECH-MEDICAL
MEASURING HUMAN PERFORMANCE



JTECH Medical Industries, Inc
7633 S Main
Midvale, UT 84047
United States of America

(385) 695-5000
(800) 985-8324
Fax (385) 695-5001
www.jtechmedical.com

